



Definition of length on YIGAO ELEPHANT II

YIGAO Sheath Specifications





Operation Tips: Insert ELEPHANT II near to stone, create high efficient suction vortex



Two steps of placement of ELEPHANT II during RIRS


- **Get to UPJ:** Placement of the ELEPHANT II to UPJ is having the same procedure like the conventional sheath. After the ELEPHANT II is placed to UPJ, the extra 5-10 cm length will be left outside urethra.
- **Reach to calyx:** Insert the ELEPHANT II together with the ureterscope to calyces, get closer to calculi. Better suction performance is guaranteed.

Possible reason: if there is always too much extra sheath is left outside during the whole surgery:

- The calculi is close to UPJ, no need to go further into calyx;
- The sheath didn't go into calyces together with scope;
- Short length of urethra.



Operation Tips: Choose proper ureterscope, keep RESD at 0.75-0.85, to keep suction performance in best status



RESD: Ration of Endoscope and Sheath Diameter

REDS= Diameter of endoscope/Diameter of UAS

When RESD is smaller than 0.85, the suction performance would be good in bendable UAS. **The recommended RESD is to keep it no bigger than 0.85.**

World Journal of Urology (2024) 42:122
<https://doi.org/10.1007/s00345-024-04815-7>

ORIGINAL ARTICLE



The optimal ratio of endoscope-sheath diameter with negative-pressure ureteral access sheath: an in vitro research

Jingyu Shi¹ · Ting Huang¹ · Baiyang Song^{1,2} · Wanzhang Liu¹ · Yue Cheng¹ · Li Fang¹

Received: 29 July 2023 / Accepted: 16 January 2024
© The Author(s), under exclusive licence to Springer-Verlag GmbH Germany, part of Springer Nature 2024

Abstract

Purpose To maintain safe intrarenal pelvic pressure (IPP), the combination of flexible ureterscope (fURS) and traditional ureteral access sheath (T-UAS) should maintain a basic rule that is the ratio of endoscope-sheath diameter ($RES D \leq 0.75$). However, the negative-pressure ureteral access sheath (NP-UAS) may break the rule of negative pressure suction. This study aimed to examine the effect of NP-UAS on IPP and flow rate (FR) with varying RESD.

Methods In a 3D-printed renal model, flexible ureteroscopy lithotripsy (fURL) was replicated. Six sizes of fURS paired with 12Fr T-UAS and NP-UAS resulted in six distinct RESDs of 0.63, 0.78, 0.87, 0.89, 0.90, and 0.91. While the irrigation pressure (IRP) was set between 100 and 800 cmH₂O and the sucking pressure (SP) was set between 0 and 800 cmH₂O, the IPP and FR were measured in each RESD.




Results NP-UASs can reduce the IPP and increase the FR at the same RESD compared to T-UASs. The IPP decreased with increasing SP with NP-UAS. When $RES D \leq 0.78$, T-UAS and NP-UAS can maintain $IPP < 40$ cmH₂O in most circumstances. When $RES D = 0.87$, it is challenging for T-UAS to sustain $IPP < 40$ cmH₂O; however, NP-UAS can do so. When $RES D \geq 0.89$, it is difficult to maintain an $IPP < 40$ cmH₂O even with NP-UAS.

Conclusion NP-UAS can decrease IPP and increase FR compared with T-UAS. To maintain a safe IPP, it is recommended that $RES D < 0.85$ when utilizing NP-UAS.

Keywords Negative-pressure ureteral access sheath · Ratio of endoscope-sheath diameter · Flexible ureteroscopy lithotripsy · Intrarenal pelvic pressure · Sucking pressure



Operation Tips: Regulate suction and irrigation flow, keep a balanced suction performance




- To maintain a balanced IRP condition and better suction performance, it can be achieved by adjusting suction and irrigation flow. To maintain a better suction performance. The recommended negative pressure parameter is: **0.02-0.024MPa/150-200mmHg**, while the irrigation flow is kept at **150-200ml/min**
- During the RIRS, it is recommended to **adjust the negative pressure and irrigation flow dynamically to maintain a balanced IRP(Intrarenal Pressure) inside kidney.**



Operation Tips: Better suction performance by placing the ureterscope to correct postion to sheath

Suggested Usage of Tip-Flexible Access Sheath



Tip of ureterscope is better to align with tip of access sheath or withdraw into tip of access sheath to realize suction.



Physical pictures under ureterscope lens, better to show edge of inner sheath tip.