

PUSEN



Müligässli 1
CH- 8598 Bottighofen

Tel. +41 71 690 00 06
www.stojanmedical.ch

 Stojan Medical

 [stojanmedical](https://twitter.com/stojanmedical)

Single Use Suction Access Ureteroscope

PU400A

World's First Certified
Single Use Suction
Access Ureteroscope



5.1 Fr

Working Channel

| Direct In-Scope Suction

| Cutting-Edge Patented Handle Design

Ergonomically designed handle for fatigue relief

| Straight Working Channel

Adjustable access port for direct instrument insertion

5.1 Fr
Working Channel

- Reduce intrarenal pressure (IRP) and temperature (IRT), mitigating the risks of postoperative infections and other complications
- Enhance Stone-Free Rate (SFR)
- Rapid cooling reduces heat damage, minimizes complications

PU400A	
Outer Diameter	9.2 Fr
Deflection	U/D 300° min.
Field of View	120°
Working Channel	5.1 Fr
Working Length	650 mm

| Compatible with PV300 Image Processor

Light and portable

Instant photo

Full-HD

PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.

5/F, Building 1, No 33, Keji San Road,
High-tech Zone, Tangjiawan Town,
519085 Zhuhai, Guangdong, People's Republic of China
(86) 400-688-2023

Rev.202503

Specifications, dimensions, appearance, and colors are subject to actual products.
Any changes will not be notified separately.