





# ELEPHANT-II

BENDABLE URETERAL ACCESS SHEATH ( NEGATIVE PRESSURE )  
G2 ADVANCED

BEND WITH THE ENDOSCOPE | NEGATIVE PRESSURE SUCTION | SUPER SMOOTH COATING

# ELEPHANT-II





**THE WORLD - FIRST INNOVATION  
BENDABLE URETERAL ACCESS SHEATH  
( NEGATIVE PRESSURE)**



## COMMENCE SUCTION WHILE STONES ARE BEING BROKEN

The sheath is adjacent to the stones, shortening the distance of the vortex flow. Greatly improve the suction efficiency of negative pressure.



# G2 ADVANCED


## NEGATIVE PRESSURE SUCTION CONNECTOR


- | MAXIMIZE STONE-REMOVAL EFFICIENCY
- | REDUCE INTRA - RENAL PRESSURE
- | CLEAR SURGICAL FIELD OF VIEW
- | AVOID THERMAL DAMAGE
- | REDUCE COMPLICATIONS

WATER SEALING CAP

PRESSURE CONTROL  
SLIDER

THE URETEROSCOPE RETRIEVAL  
INDICATOR LINE





1 Accurate scale markers  
for precise sheath placement

2 The outer sheath and the  
15cm dilator tip are covered  
with lubricious hydrophilic  
coating for providing ease  
during insertion.

3 Built-in Spring Coil  
Effectively improve the radial  
support force of the sheath.



**ZHEJIANG YIGAO MEDICAL TECHNOLOGY CO., LTD.**

Address: 3rd Floor, Building 3, Gaoxin 9th Road, Xiaoshan Economic and Technological Development Zone,  
Xiaoshan District, Hangzhou, Zhejiang

Phone: +86 571 86022937