

# activeCut



## The prostate morcellator

Experience the new AMX morcellator for anatomic endoscopic enucleation of the prostate. Designed with surgeons and procedures in mind, it combines excellent technology with an outstanding user experience.



treatment-  
induced UI



one-click collection  
system



integrated  
aspiration



integrated  
scale

The prostate morcellator activeCut enables minimally invasive treatment of **benign prostate obstruction** like never before. Thanks to its intuitive graphical user interface, simple aspiration set-up and integrated scale, the morcellator offers **unique usability** in the operating theatre.

The device features uninterrupted treatment without the need to empty a vacuum reservoir and a **reduced risk of blockage** due to fewer plug connections. While its atraumatic blade geometry ensures patient safety.

The **Active Medical Xcellence GmbH** is an German company, developing practical, secure and intuitive solutions for urology.



+49 3641 55 475 0  
[sales@amx-med.de](mailto:sales@amx-med.de)  
[www.amx-med.de](http://www.amx-med.de)

## Technical specifications

<b>activeCut (console)</b>		<b>REF 1000000</b>	<b>1</b>
device model	table top device		
aspiration	integrated aspiration pump		
classification (acc. MDR)	IIa		
display	colored 10" LCD with touch		
electrical requirements	100-240 Vac, 50/60 Hz, 120 VA		
<b>activeDrive (handpiece)</b>		<b>REF 1000002</b>	<b>2</b>
rotation speed	500–3000 rpm		
operation mode	rotation/oscillation		
sterilisation method	steam sterilisation		
cable length	3 m		
<b>MaCu OP Kit 10 (consumables)</b>			
<b>activeBlade</b> (blade)	1 x		<b>3</b>
dimensions	Ø 15 Fr x 385 mm free length		
uses	reusable (up to 10 x)		
sterilisation method	steam sterilisation		
<b>activeCollect</b> (collection system)	10 x		<b>4</b>
capacity	300 ml (300 g)		
uses	single use, non sterile		
<b>MaCu OP Kit 60 (consumables)</b>			
<b>activeBlade</b> (blade)	6 x		
<b>activeCollect</b> (collection system)	60 x		