

PUSEN

STOJAN 
Medical Systems
GmbH

Müligässli 1
CH - 8598 Bottighofen

+41 71 690 00 06
stojanmedical.ch

 Stojan Medical

 stojanmedical

Single-Use Video Ureteroscope with Direct In-Scope Suction

PU630A

6.3 Fr

DISS™



Go Further with World's First 6.3 Fr DISS™

The world's first Direct In-Scope Suction single-use ureteroscope with a groundbreaking 6.3 Fr diameter. Engineered to minimize infection risk, streamline procedures - all in a lightweight, easy-to-use design. PU630A series redefines performance in single-use endoscopy

3.6 Fr Working Channel

By maintaining the same 3.6 Fr working channel as the PUSEN 7.5 Fr ureteroscope, the PU630A series delivers uncompromised accessory compatibility

PU630A Specifications

Outer Diameter	6.3 Fr (2.1 mm)
Deflection	U/D 275°
Field of View	120°
Working Channel	3.6 Fr

With Bullet Tip

Easy insertion
Less trauma
Improved navigation



PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.

No. 1118, Jinrui Yi Road
High-tech Zone, Tangjiawan Town
519085 Zhuhai, Guangdong, People's Republic of China
(86) 400-688-2023

Rev. 202603-EN

*Specifications, dimensions, appearance, and colors are subject to actual products.
Any changes will not be notified separately.*