

# activeCut



## The prostate morcellator

Experience the new AMX morcellator for anatomic endoscopic enucleation of the prostate. Designed with surgeons and procedures in mind, it combines excellent technology with an outstanding user experience.



treatment-  
induced UI



one-click collection  
system



integrated  
aspiration



integrated  
scale

The prostate morcellator activeCut enables minimally invasive treatment of **benign prostate obstruction** like never before. Thanks to its intuitive graphical user interface, simple aspiration set-up and integrated scale, the morcellator offers **unique usability** in the operating theatre.

The device features uninterrupted treatment without the need to empty a vacuum reservoir and a **reduced risk of blockage** due to fewer plug connections. While its atraumatic blade geometry ensures patient safety.

The **Active Medical Xcellence GmbH** is a German company, developing practical, secure and intuitive solutions for urology.

## Technical specifications

<b>activeCut (console)</b>	<b>REF 100-0000</b>	<b>1</b>
device model	table top device	
aspiration	integrated aspiration pump	
classification (acc. MDR)	Ila	
display	colored 10" LCD with touch	
electrical requirements	100-240 Vac, 50/60 Hz, 90 VA	
<b>activeDrive (handpiece)</b>	<b>REF 100-0002</b>	<b>2</b>
rotation speed	500-3000 rpm	
operation mode	rotation/oscillation	
sterilisation method	moist heat (acc. ISO 17665)	
cable length	3 m	
<b>MaCu Operation Kit 10 (consumables)</b>	<b>REF 400-0003</b>	
activeBlade (blade)	1 x	<b>3</b>
dimensions	O 15 Fr x 385 mm free length	
uses	reusable (up to 10 x)	
sterilisation method	moist heat (acc. ISO 17665)	
activeCollect (collection system)	10 x	<b>4</b>
capacity	300 g	
uses	single use, non sterile	



**AMX – Active Medical Xcellence GmbH**  
 Göschwitzer Str. 22  
 07745 Jena · Germany  
 sales@amx-med.de  
[www.amx-med.de](http://www.amx-med.de)

Exclusively distributed in Switzerland by



[www.stojanmedical.ch](http://www.stojanmedical.ch)

Please note that not all products are available in all countries. Always follow the information accompanying the product. For detailed information as well as risks and side effects, please refer to the instruction for use or contact an authorized dealer or the manufacturer directly.