

PUSEN



Mülligässli 1
CH - 8598 Bottighofen

+41 71 690 00 06
stojanmedical.ch

stojanmedical


stojanmedical

Single Use Flexible Video Ureteroscope PU3033AH

Super slim 7.5 Fr outer diameter
for increased irrigation

7.5 Fr





Direct In-Scope Suction for the removal of dust and debris for enhanced visualization lower the intrarenal pressure (IRP) and temperature (IRT), reduced risks of postoperative infections and other complications, and minimized usage of accessories

PU3033AH specifications

Outer diameter	7.5 Fr (2.5 mm)
Deflection	U/D 285°
Field of view	120°
Working channel	3.6 Fr

Compatible with PV300 image processor

- HD touch screen
- Video & image capture
- Battery powered for enhanced mobility
- Integrates with surgical video systems

PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.

No. 1118, Jinrui Yi Road
High-tech Zone, 519085 Zhuhai, Guangdong
People's Republic of China
(86) 400-688-2023

Specifications, dimensions, appearance, and colors are subject to actual products. Any changes will not be notified separately.