

PUSEN



Mülligässli 1
CH - 8598 Bottighofen

+41 71 690 00 06
stojanmedical.ch

stojanmedical

stojanmedical

Single Use Flexible Video Ureteroscope PU3033A

Super slim 7.5 Fr outer diameter
for increased irrigation

7.5 Fr



| Super slim 7.5 Fr outer diameter

Compatible with a great variety of sheath types, a reduced risk of postoperative infections and other problems, and a lower RESD*

| Optimized visualization auto adjusting dual LEDs

Advanced HD image processing featuring patented image adjusting algorithms

| Compatible with 9/11 Fr UAS



Increased irrigation flow in both small and large ureteral access sheath

Model	Outer diameter	Deflection	Field of view	Working channel	Customizable buttons
PU3033A	7.5 Fr	U/D 285°	120°	3.6 Fr	2

7.5 Fr VS 9.5 Fr

Lower RESD*

Higher irrigation


Lower complication rate

Higher stone-free rate

PU3033A 7.5 Fr

Market average 9.5 Fr

12/14 Fr sheath



Global leader in single use flexible ureteroscopes

* RESD: Ratio of Endoscope-Sheath Diameter

PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.

No.1118, Jinrui Yi Road
High-tech Zone, 519085 Zhuhai, Guangdong
People's Republic of China
(86) 400-688-2023

Specifications, dimensions, appearance, and colors are subject to actual products.
Any changes will not be notified separately.