

PUSEN



Mülligässli 1
CH - 8598 Bottighofen

+41 71 690 00 06
stojanmedical.ch

stojanmedical

stojanmedical

Single Use Suction Access Ureteroscope

PU400A

World's First Certified
Single Use Suction
Access Ureteroscope

5.1 Fr_{Working Channel}



| Direct in-scope suction

| Ergonomic handle

Ergonomically designed handle for fatigue relief

| Straight working channel

Adjustable access port for direct instrument insertion

5.1 Fr

Larger working channel

- Reduce intrarenal pressure (IRP) and temperature (IRT), mitigating the risks of postoperative infections and other complications
- Enhance stone-free rate (SFR)
- Rapid cooling reduces heat damage, minimizes complications

PU400A specifications

Outer diameter	9.2 Fr
Deflection	U/D 285°
Field of view	120°
Working channel	5.1 Fr
Working length	650 mm

| Compatible with PV300 image processor

Light and portable

Instant photo

Full-HD

PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.

No. 1118, Jinrui Yi Road

High-tech Zone, 519085 Zhuhai, Guangdong

People's Republic of China

(86) 400-688-2023

Ver. 202604-EN

Specifications, dimensions, appearance, and colors are subject to actual products.
Any changes will not be notified separately.